

```

close all

%% Data Prep
% Training set
folderPath = 'train';
imageFiles = dir(fullfile(folderPath, '*.jpg'));
numImages = length(imageFiles);

for i = 1:numImages
    imageFilename = fullfile(folderPath, imageFiles(i).name);
    image = imread(imageFilename);
    trainData(i).image = image;
end
disp(trainData)
% Specify the path to the folder containing your XML files
xmlFilesFolder = 'train\';
% List all XML files in the folder
xmlFiles = dir(fullfile(xmlFilesFolder, '*.xml'));
% Initialize arrays to store bounding box information
allBboxes = [];
allLabels = {};
% Loop through each XML file
for i = 1:numel(xmlFiles)
    % Read the XML file
    xmlFile = fullfile(xmlFilesFolder, xmlFiles(i).name);
    xmlData = xml2struct(xmlFile);
    % Extract bounding box information
    annotations = xmlData.annotation;

    % Check if there is more than one object annotation
    if iscell(annotations.object)
        for j = 1:numel(annotations.object)
            % Extract label
            label = annotations.object{j}.name.Text;
            % Extract bounding box coordinates
            xmin = str2double(annotations.object{j}.bndbox.xmin.Text);
            ymin = str2double(annotations.object{j}.bndbox.ymin.Text);
            xmax = str2double(annotations.object{j}.bndbox.xmax.Text);
            ymax = str2double(annotations.object{j}.bndbox.ymax.Text);
            % Store bounding box information
            allBboxes = [allBboxes; xmin, ymin, xmax - xmin, ymax - ymin]; %
            [x, y, width, height]
            allLabels = [allLabels; label];
        end
    else
        % Extract label
        label = annotations.object.name.Text;
        % Extract bounding box coordinates

```

```

        xmin = str2double(annotations.object.bndbox.xmin.Text);
        ymin = str2double(annotations.object.bndbox.ymin.Text);
        xmax = str2double(annotations.object.bndbox.xmax.Text);
        ymax = str2double(annotations.object.bndbox.ymax.Text);
        % Store bounding box information
        allBboxes = [allBboxes; xmin, ymin, xmax - xmin, ymax - ymin]; % [x, y,
width, height]
        allLabels = [allLabels; label];
    end

end

% Get a list of all files in the folder
files = dir(fullfile(folderPath, '*.xml')); % Modify the extension based on
your image format
imageFilesNames = {files.name};
uniqueImages = unique(imageFilesNames); % Assuming imageFileNames is a cell
array of image file names
numImages = numel(uniqueImages);
boxCellArray = cell(numImages, 1);
labelCellArray = cell(numImages, 1);
for i = 1:numImages
    % Find indices corresponding to the current image
    indices = strcmp(imageFilesNames, uniqueImages{i});

    % Extract bounding box data and labels for the current image
    currentBboxes = allBboxes(indices, :);
    currentLabels = allLabels(indices);

    % Store data in cell arrays
    boxCellArray{i} = currentBboxes;
    labelCellArray{i} = currentLabels;
end

% Create a table with the required columns
boxLabelsTable = table(boxCellArray, labelCellArray, 'VariableNames',
{'BoundingBox', 'Label'});
% Create a boxLabelDatastore
boxLabelIDS = boxLabelDatastore(boxLabelsTable);
% Initialize imageDatastore
imds = imageDatastore(fullfile(folderPath, {imageFiles.name}));
% Create a boxLabelDatastore
boxLabelIDS = boxLabelDatastore(boxLabelsTable);
trainingData = (imds, boxLabelIDS);
disp(trainingData)
%% Model Arcitecture
net = resnet50;
imageSize = [416 416 3];
numClasses = 3;
anchorBoxes = [16, 16; 24, 24; 32, 32];
featureLayer = 'activation_49_relu';

```

```

lgraph = fasterRCNNLayers(imageSize,numClasses,anchorBoxes,net,featureLayer);
%% Training
% Step 7: Train the Faster R-CNN model
checkpointDir = 'C:\Users\mpetr\Downloads\BCC.v4-416x416_aug.voc\CheckPoint\';
% Create the checkpoint directory if it doesn't exist
if ~exist(checkpointDir, 'dir')
    mkdir(checkpointDir);
end
options = trainingOptions('sgdm', ...
    'MiniBatchSize', 1,...
    'MaxEpochs', 1, ...
    'InitialLearnRate', 1e-4, ...
    'LearnRateSchedule', 'piecewise', ...
    'LearnRateDropFactor', 0.1, ...
    'LearnRateDropPeriod', 3,...
    'CheckpointPath', 'CheckPoint','Verbose',true);
detector = trainFasterRCNNObjectDetector(trainingData, lgraph, options);
checkpointFile = 'detector_checkpoint.mat';
save(fullfile(checkpointDir, checkpointFile), 'detector');
trainedDetector = load(fullfile(checkpointDir, checkpointFile));
detectors = trainedDetector.detector;
% Step 8: Test the model on new images
% Use the trained model to perform object detection on new images
testImage =
imread('test\BloodImage_00359_jpg.rf.03331bc3903822adbe982fe9e7cd061b.jpg');
resizedTestImage = imresize(testImage, imageSize(1:2));
[bbox, scores, labels] = detect(detectors, resizedTestImage);

disp(imageSize)
% Step 9: Visualize the results
% Display the original image with bounding boxes around detected objects
detectedImg = insertObjectAnnotation(testImage, 'Rectangle', bbox,
LineWidth=5);
disp(bbox)
disp(labels)
imshow(detectedImg)
%% Step 10
% Model Evaluation
groundTruthData = [];
for i = 1:length(boxLabelDS)
    groundTruthData{i} = boxLabelDS.read(i);
end
truePositives = [];
falsePositives = [];
falseNegatives = [];
ious = [];
for i = 1:length(testImages)
    % Load test image
    testImage = imread(testImages{i});

```

```

    % Apply Faster R-CNN model to detect objects
    [bboxes, scores, labels] = detect(fasterRCNNModel, testImage);
    % Match predicted bounding boxes with ground truth bounding boxes
    groundTruthBoxes = groundTruthData{i}.boundingBoxes;
    matchedBoxes = matchBoundingBoxes(bboxes, groundTruthBoxes);
    % Calculate TP, FP, FN, and IoU
    tp = sum(matchedBoxes(:, 3));
    fp = sum(matchedBoxes(:, 1)) - tp;
    fn = sum(groundTruthBoxes(:, 2)) - tp;
    matchedIous = matchedBoxes(:, 4);
    truePositives = [truePositives; tp];
    falsePositives = [falsePositives; fp];
    falseNegatives = [falseNegatives; fn];
    ious = [ious; matchedIous];
end
% Calculate precision, recall, F1 score, and mean IoU
precision = truePositives / (truePositives + falsePositives);
recall = truePositives / (truePositives + falseNegatives);
f1Score = 2 * precision * recall / (precision + recall);
meanIoU = mean(ious);
% Calculate mean average precision (mAP)
map = calculateMAP(groundTruthData, testImages, bboxes, labels);
meanAP = mean(map);
% Display evaluation results
disp(['Precision: ', num2str(precision)]);
disp(['Recall: ', num2str(recall)]);
disp(['F1 Score: ', num2str(f1Score)]);
disp(['Mean IoU: ', num2str(meanIoU)]);
disp(['Mean AP: ', num2str(meanAP)]);
%% Step 11: Results
% Get Ratio of WBC and RBC
% Extract the first two labels
firstTwoLabels = labels(2:3);
% Count the number of times each label appears
uniqueLabels = unique(firstTwoLabels);
counts = zeros(length(uniqueLabels), 1);
for i = 1:length(firstTwoLabels)
    for j = 1:length(uniqueLabels)
        if strcmp(firstTwoLabels(i), uniqueLabels(j))
            counts(j) = counts(j) + 1;
        end
    end
end
end
% Calculate the ratio between the first two labels
ratio = counts(1) / counts(3);
% Print the results
if ratio < 6 %Values minimized for available data
    disp("Chance of AML: LOW");
elseif ratio < 10

```

```
    disp("Chance of AML: MEDIUM");  
else  
    disp("Chance of AML: HIGH");  
end
```